

Neuroglancer Mac Shortcut Key Reference

Initial State Commands			
Action	Key Shortcut	Single-User	Multi-User
Create annotation layer to activate drawing tools	ctrl - plus(+)	y	y
Volume tool (v). This gates all annotation functions	v	“Draw mode” = green “Edit mode” = red “Close session” = no color	“Draw mode” = green “Edit mode” = red “Close session” = no color
Access “draw mode” via sidebar	Sidebar “v” dropdown	y	y
Access “draw mode” via key-stroke	ctrl-z	y	Y (unstable between edit/draw mode)
Access “edit mode” via sidebar	Sidebar “v” dropdown	y	Y (unstable between edit/draw mode)
Access “edit mode” via key-stroke	ctrl-x	y	y(unstable)
Access “view mode” via key-stroke to pan across section	ctrl-x	Y - in “view mode” only	“view mode” is indicated but panning not functional

Cell Marking Commands			
Action	Key Shortcut	Single-User	Multi-User
Move cell marker	Single-finger click and drag annotation	Y - in “edit mode”	Y - in “edit mode”
Delete a cell marker	Select on sidebar and click trash can icon	y	y

Neuroglancer Mac Shortcut Key Reference (con't)

Polygon Creation Commands			
Action	Key Shortcut	Single-User	Multi-User
Draw a polygon	Single-finger click on trackpad	y	y
Close polygon vertex in “draw mode”	double-finger click on trackpad	y	y
Add a vertex in “edit mode”	ctrl + option + single-finger click	y	y
Delete a vertex in “edit mode”	ctrl + option + two-finger click	y	y
Pin a polygon (for cloning) MUST BE in “edit mode”	Red thumbtack in blue banner selection	y	y
Clone a polygon after pinning and MUST BE in “draw mode”	ctrl-c	y	y
Shift a polygon	shift + single-finger click and drag anywhere on polygon	Y - only in “edit mode”	Y - in “edit mode” or “draw mode”
Rotate polygon left	r	Y - only in “edit mode”	y
Rotate polygon right r	e	Y - only in “edit mode”	Y - only in “edit mode”